

Pütter® Flex

Bi-elastische Kurzzugbinden
Bi-elastic short-stretch bandages
Bandes bi-extensibles à allongement court
Bi-elastische zwachtels met korte rek
Bende bielastiche a corta trazione
Vendas bielásticas de corta elasticidad
Ligaduras bielásticas de curta tração
Επίδεσμοι μικρής τάνυσης με οριζόντια και κάθετη ελαστικότητα
Opaski bi-elastyczne o niskiej rozciągłości
Бинты короткой растяжимости, способные растягиваться в обоих направлениях
Bandaje bi-elastice cu elasticitate redusă
Bielastiska kortsträcksbindor
Bielastiset vähäelastiset siteet
雙彈性短拉力繃帶
أربطة مزدوجة المرونة قصيرة المدى

1

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EN Instructions for use

General description of the device

PütterFlex* is a brown-coloured, bi-elastic short-stretch bandage with textile elasticity and an extensibility of approximately 90 % lengthwise and approximately 40% widthwise for very high compression.

Intended purpose

PütterFlex can be used for long-term compression therapy:

- treatment of chronic venous insufficiency (CVI) with and without venous leg ulcers (VLU)
- thrombosis / thrombophlebitis
- chronic oedemata of other origins (ankle brachial pressure index – ABPI 0.8–1.3)
- lymph oedemata

The product can be used on humans by healthcare professionals in clinical and home environment in combination with padding bandages and inner lining material. PütterFlex is reusable.

Properties

Extensibility: approx. 90% lengthwise, approx. 40% widthwise.

Washable: Wash the bandage in a regular wash cycle at a temperature of 95°C.

Indication

The product can be used for the treatment of varicose veins, thrombosis / thrombophlebitis, chronic venous insufficiency (CVI), primary and secondary prevention of leg ulcer venous/venous insufficiency and oedemata of other origins (ankle brachial pressure index - ABPI 0.8 - 1.3).

Intended users

PütterFlex is designed to be used by healthcare professionals who have received appropriate training. PütterFlex can be used for the treatment of: varicose veins, thrombosis / thrombophlebitis, post-thrombotic syndrome, thromboprophylaxis, chronic venous insufficiency (CVI), venous insufficiency, angiodysplasia and other oedemata.

Composition

100 % cotton

Contraindications

PütterFlex must not be used in case of:

- advanced peripheral arterial occlusive disease (PAOD)
- decompensated heart failure (NYHA II +IV)
- phlegmasia coerulea dolens
- in case of ankle-brachial pressure index (ABPI) > 1,3 and \leq 0,5
- septic phlebitis

Precautions

Special precautionary measures are required in the event of:

- pronounced sensibility disorders of the extremities (e.g. in the case of advanced diabetic peripheral polyneuropathy)
- pronounced weeping dermatoses
- incompatibility with compression bandage material
- primary chronic polyarthritis

The application technique may have to be adapted to the patient's anatomical and pathophysiological characteristics (e.g. in regard to stretch, padding, etc.). If the foot pulse is not palpable, a Doppler ultrasound should be carried out to determine the ABPI prior to starting therapy. If an inelastic system is created by close to full stretch application, high ankle pressures should be closely monitored by a physician.

Application (see illustrations)

Modified bandaging technique according to Pütter, if not otherwise prescribed by the doctor. (1)

Further remarks

PütterFlex is indicated for use by healthcare professionals who have received appropriate training. Before applying PütterFlex, clean the wound and cover it with a suitable wound dressing. If required, excess padding bandage material may be used for additional padding of bony prominences and tissue recesses (e.g. shinbone, Bisgaard's backdrop). If additional fixation is required, use tape but no bandage clips. After application, the compression bandage must be checked (e.g. circulation in the patient's toes, pain reported by the patient) in order to be certain that the pressure of the bandage is not too high. If at any time during therapy with PütterFlex the patient develops severe pain, numbness, tingling or noticeable discolouration of the toes, a healthcare professional must be contacted immediately and the bandage must be removed. During therapy with PütterFlex, the patient should be encouraged to move.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with regular hospital waste.

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